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**THINH PHAT SCALES JOINT STOCK COMPANY**

Address : 57 Street D1, Ward 25, Binh Thanh District

HoChi Minh City Viet Nam

Tel : (08) 62.999.111 , (08) 62.888.666 , (08) 3512.7131(12 line)

Fax : (08) 3512.7699 , Order (08) 62.999.666

Website : [www.canthinhphat.com](http://www.canthinhphat.com) - [www.canvina.com](http://www.canvina.com)

Email : [hcm@canvina.com](mailto:hcm@canvina.com) - [vo@canvina.com](mailto:vo@canvina.com)

## I. CAUTIONS AND PREPARATION

### CAUTION

- 1.1 Please turn off the charge mode if do not use the hydrogens battery of nickel or Lead — aid battery.
- 1.2 Avoid installing the balance in direct sunlight-this may cause poor performance or malfunctions.
- 1.3 Do not mix different type of batteries. Replace all the Batteries at the same time.
- 1.4 If the balance is not to be used for a long period of time, remove all batteries from the battery compartment to avoid leakage which may cause damage to the instrument.
- 1.5 Avoid overloading or dropping material onto the platform as that could damage the balance.
- 1.6 Do not spill liquids on the balance as it is not water-resistant. Liquids may damage the case and if it gets inside the balance it may cause damage to the electronics.
- 1.7 Material that has a static electric charge could influence the weighing. Discharge the static electricity of the samples, if possible. Another solution to the problem is to wipe both sides of the pan and the top of the case with an antistatic agent.

### PREPARATION

- 1.1 Place your balance on a firm and level table and then level the balance by adjusting the feet to center the bubble in the spirit Level.
- 1.2 Avoid placing the balance in direct strong sunlight or draught that may cause poor performance or malfunctions.
- 1.3 Use the individual socket to avoid the interfere of other wirings.
- 1.4 Turn on the balance while there is no load on the pan.
- 1.5 Place the weight on the center of the pan for accurate weighing. Make sure the weight do not over the Max capacity as it could damage the load cell inside.
- 1.6 We suggest to warm — up the balance for 5 minutes for accurate weighing.
- 1.7 Calibration may be required before weighing or when the balance is initially installed


or moved from a location.

1.8 Please replace the battery when the low battery  mark comes.

## 2.0 SPECIFICATIONS


Model	VIBRA TPS3	VIBRA TPS6	VIBRA TPS15	VIBRA TPS30
Certificate Number	1047-2014	1048-2014	1049-2014	1050-2014
Maximum Capacity	3 kg	6 kg	15 kg	30 kg
Division (d)	0.1 g	0.2g	0.5 g	1 g
Readability (e)	0.5 g	1g	2 g	5 g
Min	10g	20g	40g	100g
Display Resolution	1: 30000	1: 30000	1: 30000	1: 30000
Non-linearity	3~5 seconds			
Display	LCD 6 number Model: TPS Seri C_DOUBLE_HT1621LCD DB05053A Revision: 1.1 produce by TPS Corporation			
Loadcell	Model: LAB-B (B is capacity of loadcell) produce by TPS Corporation			
Unit	kg, g, lb and oz			
Overload	Balance Scales will automatically protect overload exceed 9d when compared with the largest balance is installed			
Safe overload	150 %			
Power source	Pin -DC6V 4Ah Adaptor -input: AC 110-240V output: 12V/500mA			
Temperature	preservation: -10°C~+50 °C; Work : 0°C~+40 °C			
Moisture	preservation: 5%R.H~90%R.H; Work: 10%R.H~80%R.H			
Dimensions (LxDxH)	(264 x 315 x 110) mm			
Pan-size (LxD)	(260 x 218) mm			
Net Weight	3.6/4.4 kg			

Dimension(mm):

### 3.0 DISPLAY AND KEY DESCRIPTION

LED display

## VIBRA SHINKO

Indication sign: stable zero Tare low battery

	UNIT shift key in weighing mode Kg, g, lb and oz.		MODE Shift weighing and counting function on weighing mode In setting mode, it could shift the selection
	ZERO Zero in weighing mode Add key in counting mode.		TARE Tare key in weighing mode Reduce key in counting mode

#### 4.0 OPERATION INSTRUCTION

4.1 Power on Turn on the power, and the scale auto — check from 0 to 9 sequentially. Then it will show voltage of  for 2 second and go to weighing status.

4.2 Zeroing If the value of weighing is not zero (0 is not appears) before weighing, press [ZERO], 0 will come. Note: the zero range is 4% of max capacity.

4.3 Units Selection: Select the current unit in the weighing mode by press [UNIT] key.

4.4 Tare Function: Put the container on pan, press [TARE] key, the display will show the tare light.

A. Move the container it will show minus weight.

B. Press [TARE], it will show zero then gone.

4.5 counting function:

a. in weighing state, press a sample weight (the weight should over 1% of the Max capacity.) Press [mode] go to the counting setting mode, then it will display the sample number.

b. Press [Zero] to add the sample number and [Tare] to reduce the number. 10 pcs will be added or reduced for one time. If you hold the key for one second, the number will change one by one.

c. confirm the sample weight the press [mode] go to counting status. Now it will display the actually weight number.

d. press [mode] to go back to weighing status.

**4.6 Communication of serial interface(optional)** The communication of serial interface contains four modes.

a. Continuous transfer: don't receive orders and transfer the data every 0.3 seconds in the weighing mode.

- b. Response transfer: the pc transfers the orders of inquiring and zeroing, the scale will response and transfer the corresponding datas.
- c. Stable transfer: don't receive the orders and transfer the data once the weighing is stable.
- d. Print mode: The mode is same as the mode 3) but the transferng format is different. It applies to printer connection. The detailed setting refers to 6.3 Communication interface setting.

## 5 .0 ALARMING FUNCTION



### 5.1 Over load alarming

When the load on the platform is over 9 D, the beep will sound and the display will show-----, you must remove the load from platform immediately.

### 5.2 Low load alarming:

If there is problem in load cell, the beep will sound and the display will Show  

### 5.3 Counting alarm:


In counting status, when the scale is in minus number it will display  and when the number is over 1,000,000, it will display 

**5.4 Low battery alarming:** When the light of low battery voltage appears, you're requested to recharge the scale over 12 hours. When the scale is in recharging, the mark light will be in flashing. The mark light will disappear after recharge have finished.


**5.5 Charge indication:** When the scale was charging, the charge light is red. It will turn green when finishing the charge.

**6 .0 SETTING INSTRUCTION:** Hold [TARE] and turn on the scale to enter the setting mode. It will show the selection of parameter. Press [Mode] to go to the selection.

**6.1 Unit setting:** This function is for the setting of unit mode. The unit can't be shown on the scale if it is turned off by this setting in weighing mode. Press  to select the units, press  to turn on or off the units. The details are as follow:


**6.2 Backlight setting:** The function is for the setting of backlight. When the display shows b1 2, press [ZERO] to select the backlight mode.


**6.3 Communication Interface setting:** Communication interface setting includes sending mode, address and BPS selection. Please refer to 7.1 communication protocols in addenda to see the details of protocols.

**Communication mode:** This function is for selecting communication mode. There are four modes selectable: continuous transfer mode, response transfer mode; stable transfer mode and printing mode. Press [ZERO] to select mode and press [UNITS] to enter the selection of indicator address. Enter the unit setting by pressing [MODE]. The details are as follow:


**6.4 Selection of indicator address:** This parameter is for selecting indicator address. There is a data for indicator address in the communication protocols. You can distinguish the different indicator by the address. The selection range is from 1- 25. Please refer to 5.1 interface communication protocols in addenda to see the details of protocols. Press [ZERO] to select address or [UNITS] to enter the selection of BPS. Press [MODE] to enter the unit setting. he details are as follows:


**6.5 Selection of BPS :** This is used to select BPS. When BPS of the indicator and the communication equipment is not same, the communication can't be used normally. Press **ZERO** to select BPS or **UNITS** to enter the selection of communication mode. Press **MODE** to enter the unit setting. The details are as follows:

## 7.0 ATTACHMENT

**7.1 The balance use RS 232 UART signal, one frame consists of 10 Bit.**

BIT1: Start bit BIT2 — BIT9: Data bit BIT10: Stop bit.

**The format of printing mode(mode 4) is as follow:**

**add: A N0.1 : Weighing times**  
**n/w: 0.00001 kg: net weight**

**The format of continue transfer mode (mode 1) is as follow:**

All data is on ASCII Code except the start bit.

Byte	Notes	Contents
1	Start of text	02H(STX)
2	Indicator Add	A~Z
3	Stable status	0 stable/1 unstable
4	Data Sign	+/-
5~10	Data Sign	6 Byte (From high to low)
11	Decimal	From right to left (0-5).
12~14	Unit Code	A~Z
15	Power	0 normal/1 low battery
16, 17	XOR Checksum	2 Byte A~Z
18	End of text	03(ETX)

XOR=28 3\*\*\*\*\* 15

Note: Checksum of high — half byte and low-half byte less or equal to 9, adding 30H, then the data transfer as number of ASCII code. Example: Checksum of high — half is 6, adding 30H, it come up to 36H, data will be transferred as 6 of ASCII code. Checksum of high — half byte and low- half byte is over 9, adding 37H, the data transfer as alphabet of ASCII code. Example: Checksum of high — half is B, adding 37H, it come up to 42H, the data will be transferred as B of ASCII code.

**7.2 Unit Conversion (The units used may be different in different models)**

1kg	=	1000g
1lb	=	453.59237g
1oz	=	28.349523125g
1HL	=	605g
1t(TOLA)(INDIA)	=	11.6638038g

**Hotline Call to our company TPS Corporation in Viet Nam**

Call	Tel
Hotline	(+84.8)62.888.666 / (+84.8) 62.999.111
Mr Cuong (24/24)	0084 915.999.111
Mr Xa	0084 974.000.333
Ms Cuong	0084 908.444.000

**THINH PHAT SCALES JONT STOCK COMPANY**

Address : 57 Street D1, Ward 25, Binh Thanh District  
HoChi Minh City Viet Nam

Tel : (08) 62.999.111 , (08) 62.888.666 , (08) 3512.7131 (12 line)

Fax : (08) 3512.7699 , Order (08) 62.999.666

Website : [www.canthinhphat.com](http://www.canthinhphat.com) - [www.canvina.com](http://www.canvina.com)

Email : [hcm@canvina.com](mailto:hcm@canvina.com) - [vo@canvina.com](mailto:vo@canvina.com)