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### CALIBRATION.

#### Using an External Calibration Weight

Calibration is required when the weighing machine is initially installed or if the balance is moved to a substantial distance from the original location or if ambient conditions change considerably. 30 minutes' warm-up of the balance is always needed before weighing and calibration.

- Turn the Balance on and when display shows “-----”, press ZERO key
- The display will show “CAL”
- Press ZERO key again, the display will show “X000.00”, “X” is the flashing digit
- Press UNIT key, the flashing digit will move to right

- Press the TARE (or PCS, depending on the model) key to increase the value of the flashing digit (“X” is the calibration weight can be set according to users’ requirements, we recommend a minimum weight of at least 50% of the Balance capacity)
- Press the ZERO key, the display will show a series of digits known as the A/D value, these numbers will continuously move, once the display shows the stable indicator (please see page 18 for help locating the stable indicator), press the UNIT key, the display will then show “X00.000” or “X000.00” depending on the model and capacity.
- Place the corresponding test weight on the center of the pan
- The display backlight lights up and returns to weighing mode
- Now the calibration is completed
- Turn off the power, then turn on the power again, place a known weight on the platform to verify if it is accurate.
- If the weight displayed does not correspond with the value of the calibration weight, please repeat above steps.