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High Precision Balance Service Manual

DS seri

1. The function will be effective after you press and release the key.
2. when the balance is in the weighting mode, Press [mode] key for 3s to enter the parameter setting, it will show the first parameter setting "P1 UNT". Always press [mode] key to enter the next main-item. Press [Print] key to enter the current sub-item which in main item, saving and enter to next sub-item setting. [unit] key is use to modify one digit parameter, or adjust the digit for multi-digit parameter.,[tare] key is add up, [zero] key is reduce.
3. Detailed description as follows:

Main-item	Sub-item	function	Way of adjust
P1 UNT (UNIT Settings)	Kg XXX	To set the unit which is effective, XXX is "ON" or "OFF". LED model can only show 4 units(different model may shows different units),so only the first four units which set "on" can be effective and shown on the display. LCD models every unit which set "on" can be effective and shown on the display.	Press [unit] to set whether the unit is "ON" or "OFF", press [Print] to enter the next setting.
P2 TRA (Communication Settings)	FC X	Communication mode: 1.continuous transfer mode 2. Stable transfer mode 3. Response transfer mode.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	Add XX	Communication Address, effective in multi-machine communication.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	b XXXX	Communication baud rate.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	X X X	The parameter can be set to: n81 check bit/o71 data bit /e71 stop bit.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	CC XXX	Turn on/off XOR check.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
P3 FUN (Other function setting)	Pr XXX	Setting the parameter of sleep, could set to OFF (it is no sleep), or you can set every 6s as a step, setting 6-60 seconds for sleep time. LCD model is for the backlight.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.

	LtN X	X is the level of brightness. This parameter is only for LED model.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	BL XX	XX is "ON"/"OFF" of backlight, This parameter is only for LCD model.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	Hd XXd	Dynamic weighing setting set "OFF" means this function is off. if it is set "ON", you can set 10 d as a step, 10 d to 70 d dynamic range can set by step. Such as: set as 20d,when weight changes in 20d range, the software will calculate a suited average value, then lock, the same time a long beep go off. If weight over this range at this time, it will be unlocked automatically, and with a short beep. ---d, it means weigh must return to 0, then the weight can be unlocked, but the weight change should not exceed 80d extent.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	MLE XX	XX can be set to "ON"/"OFF", it means it the balance alarm when the weight is negative when user use the upper limit and lower limit alarm function.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
P4 HL (Weigh check setting)	-H-L -	Upper limit and lower limit alarm, "H" means upper limit range, "L" means lower limit range, which place's line is blinking means the weight in this area the alarm beeps, such as: if left line is blinking, it means the weight over the upper limit, the alarm beeps. if the place no line, it means the function is off.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	HXXXXX	This is upper limit range parameter.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit, When finished, press [print] to enter the next parameter
	LXXXXX	This is lower limit range parameter.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change

			the value of the digit. When finished, press [print] to back to the weighting mode.
PIN_	<p>Input code, if you put the right code you can enter into the follow parameter setting. If not, it will return to main-item setting of units(P1 UNT), when you input over, please press [print] to confirm.</p> <p>The original code is [units] [units] [tare] [units]</p> <p>If input this code cannot enter, please connect the factory.</p>		
P5 FIL { Filter range }	FI X	Filter intensity parameter, the greater the X value, the greater the intensity, the date shows more stability, but the date shows more slowly.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	ZE X	Zero track parameter, when zero is easy to move, you can set the parameter be greater. this parameter greater, the zero is more stability, but the small weigh value is difficult to identify.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	ZTR X	Creep tracking parameter, when load cell's output creeps, you can set this parameter. the greater the effect better, but the small weigh value is difficult to identify.	Press [units] to rework the parameter. Press [Print] to back to the current Main-item.
P6 C-D{Capacity&Division}	CXXXX	XXXX is weigh capacity, when weigh is more than +9 d, it will alarm over load.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit, When finished, press [print] to enter the next parameter.
	DXXXX	Division setting is the minimum resolution.	Press [units] to rework the parameter. Press [Print] to back to the current Main-item.
P7 LX {multi-range calibration }	<p>You can enter the multi-range calibration When X is 1(ON).when X is 0, means OFF. Note.: If you don't want to enter the multi-range calibration, please set X off, if not it will effect the line performance of the scale.</p>		Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	LINEX	The state of multi-range calibration, X is 1,2,3, it means which step worked.	Auto to enter the next parameter.
	XXXXX	<p>Value of weight, when you finished setting, please put the calibration weight.</p> <p>Remark: the second line of weight value is greater than the fist one, and the third line of weight value is larger than the second.</p>	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit,. When finished, press [print] to enter the next parameter.
P8 CAL	EXXXX	XXXX is the value of calibration weight.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit,. When finished, press [print] to enter the

			next parameter.
	XXXXX	When it is Zero AD value, it means there are not put anything on the pan.	Press [print] to calibration, press [units] enter to setting the AD enlarged multiple
	XXXXX	AD value, please put the weight on, and then waits for stable light turn on.	After the display stable, press [units] or [print] to calibration.

MAIN FEATURES AND SPECIFICATION

CAUTION

- 1.1 Please turn off the charge mode if do not use the hydrogens battery of nickel or Lead — aid battery.
- 1.2 Avoid installing the balance in direct sunlight-this may cause poor performance or malfunctions.
- 1.3 Do not mix different type of batteries. Replace all the Batteries at the same time.
- 1.4 If the balance is not to be used for a long period of time, remove all batteries from the battery

compartment to avoid leakage which may cause damage to the instrument.

1.5 Avoid overloading or dropping material onto the platform as that could damage the balance.

1.6 Do not spill liquids on the balance as it is not water-resistant. Liquids may damage the case and if it gets inside the balance it may cause damage to the electronics.

1.7 Material that has a static electric charge could influence the weighing. Discharge the static electricity of the samples, if possible. Another solution to the problem is to wipe both sides of the pan and the top of the case with an antistatic agent.

PREPARATION

1. Place your balance on a firm and level table and then level the balance by adjusting the feet to center the bubble in the spirit Level.

2. Avoid placing the balance in direct strong sunlight or draught that may cause poor performance or malfunctions.

3. Use the individual socket to avoid the interfere of other wirings.

4. Turn on the balance while there is no load on the pan.

5. Place the weight on the center of the pan for accurate weighing. Make sure the weight do not over the Max capacity as it could damage the load cell inside.

6. We suggest to warm — up the balance for 5 minutes for accurate weighing.

7. Calibration may be required before weighing or when the balance is initially installed or moved from a location.

8. Please replace the battery when the low battery  mark comes.

Introduction : This is a multi-function weighing balance that offer flexibility together with good features!

Main features:

A: Automatic zero point tracking function.

B: With multi-rang calibration.

C: Tare function.

D: Many units can be selection.

E: Simple counting function and percentage function are applicable. F: Fast, stable and accurate display reading with simple operation.

PRECAUTIONS:

1. Full charge the battery after unpacking the scale.

Recharge the battery: When battery symbol appears on the LCD display, charge the battery with AC adapter, the balance of charge will light up in red, when it becomes green means charge completed. (It takes about 8 hours to full charge the battery.)

2. Install the equipment on a level and stable surface.

3. Do not install the equipment near the air conditioning or a vibrating machine.

4. Install the equipment in an environment with steady temperature(0°C ~40°C), prevent from rapid temperature changes.

5. Independent AC adapter for this equipment is recommended, check the

6. Warm up the equipment for 15 minutes before use.
7. Avoid overloading or dropping material onto the platform.

SPECIFICATIONS for size 120mm

Model	DS223	DS323	DS302	DS602	DS1202	DS1502
Maximum Capacity	220 g	320 g	300 g	600 g	1200 g	1500 g
Division (d)	0.001 g	0.001 g	0.01 g	0.01 g	0.01 g	0.01 g
Readability (e)	e = 10 d					
Min	Min = 20 e					
Display Resolution	1/60000					
Response	3~5 seconds					
Display	LCD display with 6 digits / white backlight Model: TPS Seri C_DOUBLE_HT1621LCD DB05053A Revision: 1.1 produce by TPS Corporation					
Loadcell	Model: LAB-B (B is capacity of loadcell) produce by TPS					
Unit	g, ct, lb, oz, tlt (can be set by soft)					
Overload	150 %					
Appearance	ABS, with stainless steel pan, breeze shield					
Calibration weight	3g – 1500g					
Power source	Battery –DC6V 1.3Ah Adaptor –input: AC 110-240V output: 12V/500mA					
Temperature	preservation: -10°C~+50 °C; Work: 0°C~+40 °C					
Pan-size	Ø 100 mm			Ø 120 mm		
Dimensions	(182 x 283 x 305) mm			(182 x 283 x 70) mm		
Moisture	preservation 5%R.H~90%R.H; Work: 10%R.H~80%R.H					
Net Weight	4,0 kg			2,2 kg		

Connect the balance with PC or printer, this function is effective. Under weighing mode, press the **[PRINT]** key to print out percent weight value.


SPECIFICATIONS for size (160 x 160) mm

Model	DS2002	DS3002	DS6002	DS3001	DS6001	DS7501
Maximum Capacity	2 kg	3 kg	6 kg	3 kg	6 kg	7.5 kg
Division (d)	0.01 g	0.01 g	0.05 g	0.1 g	0.1 g	0.1 g
Readability (e)	e = 10 d					
Min	Min = 20 e					
Display Resolution	1/60000					
Response	3~5 seconds					
Display	LCD display with 6 digits / white backlight Model: TPS Seri C_DOUBLE_HT1621LCD DB05053A Revision: 1.1 produce by TPS Corporation					
Loadcell	Model: LAB-B (B is capacity of loadcell) produce by TPS					

	Corporation
Unit	g, ct, lb, oz, tlt (can be set by soft)
Overload	150 %
Appearance	ABS, with stainless steel pan, breeze shield
Calibration weight	3g – 1500g
Power source	Battery –DC6V 1.3Ah Adaptor –input: AC 110-240V output: 12V/500mA
Temperature	preservation: -10°C~+50 °C; Work: 0°C~+40 °C
Pan-size	(160 x 160) mm
Dimensions	(182 x 283 x 70) mm
Moisture	preservation 5%R.H~90%R.H; Work: 10%R.H~80%R.H
Net Weight	2.3 kg


Connect the balance with PC or printer, this function is effective.

Under weighing mode, press the **[PRINT]** key to print out percent weight value.

DISPLAY AND KEY DESCRIPTION

LCD display

Indication sign: stable zero Tare low battery

Start to use: Accessories:

1. adapter
2. balance body
3. plastic pan
4. SS pan
5. wind shield

Turn on the indicator

- Press the ON/OFF] key.
- Wait till the screen displays “0.00”. If it doesn’t shows “0.00”, press [ZERO] key to enter weighing mode.

KEYPAD AND PANEL INTRODUCTION

	Lock the display weight value		Change Unit g, ct, lb, oz, tlt
	Change Weight/PCS/% functions		TARE
	Print out the data		Reset to zero

FUNCTIONS AND OPERATION

1. Simple weighing

1. Press the on/off] key or Turn on the power switch on the right side of balance.

2. The balance power on with beep and backlight flash.
3. Wait till the screen displays “0.00g”. Press **[ZERO]** key to enter weighing mode if the screen doesn’t display “0.00g”.
4. Put the object on weighing pan then read the weight within 2~3 seconds.

2. Tare function:

- Put the container (object) on weighing pan, press **[TARE]** to deduct the container(object) weight, the screen displays “0.00g”.
- Put more objects on the weighing pan, the screen displays weight of new objects excluding weight of container.

3. Units shift and setting

Under weighing mode, press **[UNITS]** key to change to the unit you need. Open more units or close any units, please contact your local agent.

4. Simple counting setting

1. Under weighing mode, press **[PCS]** choose the simple counting function. Put on the weights, if no unit weight setting before, screen will display “0 pcs”.
2. Press **[UNIT]** key, enter quantity setting, the screen displays “N-XXX” (“XXX” stands for sample quantity.)
3. Press **[UNIT]** key, move the digit.
4. Press **[TARE]** key, add up the percent digits.
5. Press **[ZERO]** key, reduce the percent digits.
6. Press **[PCS]** key, save the above setting, you can begin to use this function.

5. Simple percent counting

In counting mode, press **[MODE]** to change % mode or weighting mode.

- Under weighing mode, press [MODE], choose the % function. Put on the weights, screen will display “ S-100 %”.
- Press [MODE] key, save the above setting, you can begin to use this function.

** LED model do not have % mode, so only change to weighing mode.

6. Hold function

1. Press once to lock the percent weight value.
2. Press again to cancel lock status.

** This function usually set to close, if you want to use this function, please contact your local agent.

7. Printing function

connect the balance with PC or printer, this function is effective.

Under weighing mode, press the [PRINT] key to print out percent weight value.

8. Simple calibration

	<p>Press hold key 2 seconds to enter simple calibration</p>	<p style="text-align: center;">CAL</p> <p>Press [UNIT] to move the digits Press [TARE] key, add up the percent digits Press [ZERO] key, reduce the percent digits Put on the right weights (the same as the display weight).</p>
	<p>Press key</p>	<p>Wait the stable symbol appears for 2 seconds the buzzer will send out a short beep, that means simple, calibration successfully.</p>

ERROR MESSAGE AND TROUBLESHOOTING

Error message

Message	Problem
ADO---	Over A/D resolution range
ADL---	Low A/D resolution range
OVR---	Over load (max: capacity +9e)
Battery symbol	Low battery alarm
Lo-bAt	Low battery, Need charge

Trouble shooting

Where	Possible Error	Trouble shooting.
-------	----------------	-------------------

	message	
Power on	ADO---	Check and remove the object from weighing pan or send back for
Power on	ADL---	send back for repair.
Power on	Battery symbol	Charge the battery.
Normal weighing mode	OVR---	Check if weighing object over the capacity +9e.
Normal weighing mode	Lo-bAt	Charge the battery.

High Precision Balance Service Manual

DS seri TPS

1. The function will be effective after you press and release the key.
2. when the balance is in the weighting mode, Press [mode] key for 3s to enter the parameter setting, it will show the first parameter setting "P1 UNT". Always press [mode] key to enter the next main-item. Press [Print] key to enter the current sub-item which in main item, saving and enter to next sub-item setting. [unit] key is use to modify one digit parameter, or adjust the digit for multi-digit parameter.,[tare] key is add up, [zero] key is reduce.
3. Detailed description as follows:

Main-item	Sub-item	function	Way of adjust
P1 UNT (UNIT Settings)	Kg XXX	To set the unit which is effective, XXX is "ON" or "OFF". LED model can only show 4 units(different model may shows different units),so only the first four units which set "on" can be effective and shown on the display. LCD models every unit which set "on" can be effective and shown on the display.	Press [unit] to set whether the unit is "ON" or "OFF", press [Print] to enter the next setting.
P2 TRA	FC X	Communication mode: 1.continuous transfer mode	Press [units] to rework the parameter. Press

(Communication Settings)		2. Stable transfer mode 3. Response transfer mode.	[Print] to enter next sub-item to set the detail parameter.
	Add XX	Communication Address, effective in multi-machine communication.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	b XXXX	Communication baud rate.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	X X X	The parameter can be set to: n81 check bit/o71 data bit /e71 stop bit.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	CC XXX	Turn on/off XOR check.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
P3 FUN (Other function setting)	Pr XXX	Setting the parameter of sleep, could set to OFF (it is no sleep), or you can set every 6s as a step, setting 6-60 seconds for sleep time. LCD model is for the backlight.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	LtN X	X is the level of brightness. This parameter is only for LED model.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	BL XX	XX is "ON"/"OFF" of backlight, This parameter is only for LCD model.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.

	Hd XXd	Dynamic weighing setting set "OFF" means this function is off. if it is set "ON", you can set 10 d as a step, 10 d to 70 d dynamic range can set by step. Such as: set as 20d,when weight changes in 20d range, the software will calculate a suited average value, then lock, the same time a long beep go off. If weight over this range at this time, it will be unlocked automatically, and with a short beep. ----d, it means weigh must return to 0, then the weight can be unlocked, but the weight change should not exceed 80d extent.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	MLE XX	XX can be set to "ON"/"OFF", it means it the balance alarm when the weight is negative when user use the upper limit and lower limit alarm function.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
P4 HL (Weigh check setting)	-H-L -	Upper limit and lower limit alarm, "H" means upper limit range, "L" means lower limit range, which place's line is blinking means the weight in this area the alarm beeps, such as: if left line is blinking, it means the weight over the upper limit, the alarm beeps. if the place no line, it means the function is off.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	HXXXXX	This is upper limit range parameter.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit, When finished, press [print] to enter

			the next parameter
	LXXXXX	This is lower limit range parameter.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit. When finished, press [print] to back to the weighting mode.
PIN_	Input code, if you put the right code you can enter into the follow parameter setting. If not, it will return to main-item setting of units(P1 UNT), when you input over, please press [print] to confirm. The original code is [units] [units] [tare] [units] If input this code cannot enter, please connect the factory.		
P5 FIL { Filter range }	FI X	Filter intensity parameter, the greater the X value, the greater the intensity, the date shows more stability, but the date shows more slowly.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	ZE X	Zero track parameter, when zero is easy to move, you can set the parameter be greater. this parameter greater, the zero is more stability, but the small weigh value is difficult to identify.	Press [units] to rework the parameter. Press [Print] to enter next sub-item to set the detail parameter.
	ZTR X	Creep tracking parameter, when load cell's output creeps, you can set this parameter. the greater the effect better, but the small weigh value is difficult to identify.	Press [units] to rework the parameter. Press [Print] to back to the current Main-item.
P6 C-D{Capacity&Division}	CXXXX	XXXX is weigh capacity, when weigh is more than +9 d, it will alarm over load.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit, When finished, press [print] to enter the next parameter.
	DXXXX	Division setting is the minimum resolution.	Press [units] to rework the parameter. Press [Print] to back to the current Main-item.
P7 LX {multi-range calibration }	You can enter the multi-range calibration When X is 1(ON).when X is 0, means OFF. Note.: If you don't want to enter the multi-range		Press [units] to rework the parameter. Press [Print] to enter next

	calibration, please set X off, if not it will effect the line performance of the scale.		sub-item to set the detail parameter.
	LINEX	The state of multi-range calibration, X is 1,2,3, it means which step worked.	Auto to enter the next parameter.
	XXXXX	Value of weight, when you finished setting, please put the calibration weight. Remark: the second line of weight value is greater than the fist one, and the third line of weight value is larger than the second.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit,. When finished, press [print] to enter the next parameter.
P8 CAL	EXXXX	XXXX is the value of calibration weight.	Press [units] to move the location of blinking digit, press [tare] and [zero] to change the value of the digit,. When finished, press [print] to enter the next parameter.
	XXXXX	When it is Zero AD value, it means there are not put anything on the pan.	Press [print] to calibration, press [units] enter to setting the AD enlarged multiple
	XXXXX	AD value, please put the weight on, and then waits for stable light turn on.	After the display stable, press [units] or [print] to calibration.
	ADG X	Enlarge multiple parameters; X is 1 or 2, when full-load, the scale shows AD0---, please set the amplification parameter to the smaller one.	Press [units] to rework parameter, press [print] return to zero AD.
P9 EC	E1 Rul (1s later) XXXX	Wait for 1s, please note down the value of XXXX. When the scale needs calibration, and there is no weight to use, please enter the value, then the scale can be use again.	If you want to change this parameter, first, press [Unit],at this time one of the digit will blinking, then press [tare] and [zero] to change the value of the digit. Press [Print] to enter the next parameter setting.
	EX XXX(1s later) XXXX	Wait for 1s, the value is multi-range calibration's value.4 value in total, if you have set multi-range calibration, the 4 value is effect. When the scale needs multi-range calibration, and there is no weight to use, please enter the value, then the scale can	If you want to change this parameter, first, press [Unit], at this time one of the digit will blinking, then press [tare] and [zero] to change the value of the digit. Press [Print] to enter the next parameter setting.

		be use again.	
--	--	---------------	--

4. Easy calibration

Simple calibration, in order to use more convenience, this product has simple calibration function. In weighing mode, take the weights down, and press [mode] for 3 seconds, the display will show -CAL-, then show the weight, at this time to amend the value to calibration weight, then put the calibration weight, after 3 seconds, press [print] or [mode] to calibration.

Attentions: when calibration, if the result of easy calibration compares with internal calibration are differ 100%, the easy calibration is not pass. Such as: if put 10kg weight, the internal calibration is 10kg, the internal calibration had passed, in easy calibration ,value setting is over 20kg or less then 5kg, the easy calibration will not pass.

5. Easy counting function

Easy counting function, if the unit of pcs is "on", and the PCS is in the first four units, it means the easy counting function is on.

In weighing mode, press [mode] to change the units to PCS, if you have not set the single value before, the display will show "0", at the same time put the weight on, and press [mode], the display will show " N-XXX", XXX is the sampling quantity of weight, press [units]to amend location of blinking, press [tare] and [zero] to amend location of blinking, when amend is finished, please press [mode] to sampling, at the same time the scale will counting.

In counting mode, press [mode] to change % mode or weighting mode.LED model do not have % mode, so only change to weighting mode.